

Product Datasheet

Influenza A virus(H7N9) Hemagglutinin antibody (Biotin) (orb243408)

Description	Influenza A virus(H7N9) Hemagglutinin antibody (Biotin)
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Influenza A virus Hemagglutinin protein (19-524AA)
Target	HA
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	R4NN21
Expiration Date	12 months from date of receipt.