

Product Datasheet

Lysyl oxidase homolog 2 antibody (FITC) (orb243512)

Description	Lysyl oxidase homolog 2 antibody (FITC)
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	FITC
Immunogen	Recombinant Mouse Lysyl oxidase homolog 2 protein (26-776AA)
Target	Loxl2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P58022
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.