

Product Datasheet

Protein S100A8 antibody (HRP) (orb243529)

Description	Protein S100A8 antibody (HRP)
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Mouse Protein S100-A8 protein (2-89AA)
Target	S100a8
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P27005
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.