

Product Datasheet

RANKL antibody (FITC) (orb243552)

Description	RANKL antibody (FITC)
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	FITC
Immunogen	Recombinant Mouse Tumor necrosis factor ligand superfamily member 11 protein (70-316AA)
Target	Tnfsf11
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O35235
Expiration Date	12 months from date of receipt.