

Product Datasheet

Mycobacterium tuberculosis Immunogenic protein MPT64 antibody (Biotin) (orb243567)

Description	Rabbit polyclonal antibody against Mycobacterium tuberculosis Immunogenic protein MPT64 conjugated t...
Species/Host	Rabbit
Reactivity	Protozoa
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Mycobacterium tuberculosis Immunogenic protein MPT64 protein (24-228AA)
Target	mpt64
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P9WIN9
Expiration Date	12 months from date of receipt.