

## Product Datasheet

# Mycobacterium tuberculosis Immunogenic protein MPT64 antibody (HRP) (orb243569)

<b>Description</b>	Rabbit polyclonal antibody against Mycobacterium tuberculosis Immunogenic protein MPT64 conjugated t...
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Protozoa
<b>Conjugation</b>	HRP
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Recombinant Mycobacterium tuberculosis Immunogenic protein MPT64 protein (24-228AA)
<b>Target</b>	mpt64
<b>Preservatives</b>	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<b>P9WIN9</b>
<b>Expiration Date</b>	12 months from date of receipt.