


Product Datasheet

milk beta-Lactoglobulin antibody (orb243816)

Description	milk beta-Lactoglobulin antibody
Species/Host	Rabbit
Reactivity	Bovine
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Bovine milk Beta-lactoglobulin (Native Protein)
Target	LGB
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P02754
Dilution Range	WB:1:1000-5000
Expiration Date	12 months from date of receipt.



Western blot analysis of Bovine milk Lac...