


Product Datasheet

Gliadin antibody (FITC) (orb243827)

Description	Gliadin antibody (FITC)
Species/Host	Rabbit
Reactivity	Plant
Conjugation	FITC
Immunogen	Wheat Gliadin (Native Protein)
Target	Wheat Gliadin
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.



Western blot analysis of (Lane 1): wheat...