

Product Datasheet

Corticosterone antibody (orb243880)

Description	Corticosterone antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Corticosterone-BSA conjugate
Target	Corticosterone
Preservatives	Preservative: 0.1mol/L NaCl, 15mmol/L NaN ₃ , PH7.2
Form/Appearance	liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.