



Product Datasheet

Podocalyxin antibody (orb243896)


Description	Podocalyxin antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC, IP, WB
Immunogen	Recombinant Human PODXL protein
Target	PODXL
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Uniprot ID	O00592
Dilution Range	WB:1:500-1:5000, IHC-P:1:50-1:500, IF:1:50-1:200, IP:1:200-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa(Lane 1), H...



Western blot analysis of Rat brain & hea...



Flow cytometric analysis of HeLa cells u...