

Product Datasheet

Kanamycin antibody (orb243904)

Description	Kanamycin antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Kanamycin-BSA conjugate
Target	Kanamycin
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	ELISA: 0.3-24.3 ng/ml
Expiration Date	12 months from date of receipt.