


Product Datasheet

Myoglobin antibody (orb243919)


Description	Myoglobin antibody
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human Myoglobin protein
Target	MB
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P02144
Dilution Range	WB:1:500-1:5000, IHC-P:1:50-1:500, IF:1:50-1:200
Expiration Date	12 months from date of receipt.



Western Blot analysis of Mouse heart mus...



Immunofluorescent analysis of MCF-7 cell...



Immunohistochemical of paraffin-embedded...