


Product Datasheet

Rat Camk4 protein (orb244067)

Description	Isoform 1: Calcium/calmodulin-dependent protein kinase that operates in the calcium-triggered CaMKK-...
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-tagged
Note	For research use only
Application notes	Full length of His-tag and expression region is 1-474aa
Protein Sequence	MLKVTVSPSPSSPCSSVTSSTENLVPDYWIDGSKRDPLSDFFEVESELGRGATSI VYRCKQKGTQKPYALKVLKKTVDKKIVRTEIGVLLRLSHPNIIKLKEIFETPTEISLVL ELVTGGELFDRIVEKGYYSERDAADAVKQILEAVAYLHENGIVHRDLKPENLLYAT PAPDAPLKIADFGLSKIVEHQVLMKTVCGTPGYCAPEILRGCAYGPEVDMWSVGII TYILLCGFEPFYDERGDQFMFRRILNCEYYFISPPWWDEVSLNAKDLVKKLIVLDPK KRLTTFQALQHPWVTGKAANFVHMDTAQKKLQEFNARRKLKAAVKAVVASSRL GSASSHTNIQESNKASSEQAQDQKDKTDPLENKMQAGDHEAAKAAADET MKLQSEEEVEEGVKEEEEEEEEETS RMVPQEPEDRLETDDQEMKRNSEETL KSVEEEMDPKAEAAAAGVGLGVPPQQDAILPEY
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	57.2 kDa
Uniprot ID	P13234
Expiration Date	6 months from date of receipt.



SDS-Page analysis of Rat Camk4 protein.