

Product Datasheet

Rabbit CD40LG protein (orb244085)

Description	Recombinant Rabbit CD40 ligand
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-tagged
Note	For research use only
Application notes	Expression region is 113-261aa and His-tag protein
Protein Sequence	MQKGDQDPQIAAHLISEASSKSSSVLQWAKKGYTMSNTLVTLENGKQLKVKR QGFYYIYAQVTFCSNQEPSSQAPFIASLCLKSSGGSERILLRAANARSSSKTCEQQ SIHLGGVFELQADASVFNVTASQVNHGTGFTSFGLLKL
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	20.2 kDa
Uniprot ID	G1SKP7
Expiration Date	6 months from date of receipt.