

## Product Datasheet

### Human CD44 protein (orb244086)

<b>Description</b>	Recombinant human CD44 antigen
<b>Endotoxins</b>	Not test.
<b>Preservatives</b>	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
<b>Form/Appearance</b>	Liquid or Lyophilized powder
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Tag</b>	N-terminal 6xHis-SUMO-tagged
<b>Note</b>	For research use only
<b>Application notes</b>	Extracellular domain of Isoform 4 of His-SUMO-tag and expression region is 21-606aa
<b>Protein Sequence</b>	<p>QIDLNITCRFAGVFHVEKNGRYSISRTEAADLCKAFNSTLPTMAQMEKALSIGFET            CRYGFIEGHVVIPRIHPNSICAANNTGVYILTSNTSQYDITYCFNASAPPEEDCTSVT            DLPNAFDGPITITIVNRDGTTRYVQKGEYRTNPEDIYPSNPTDDDVSSESSERSST            SGGYIFYTFSTVHPIDEDSPWITDSTDRIPTATSSNTISAGWEPNEENERDR            HLSFSGSGIDDEDFFISSTISTTPRAFDHTKQNDWTQWNPESHNPVLLQTTT            RMTDVRNGTTAYEGNWNPEAHPPLIHHEHHEEEETPHSTSTIQATPSSTTEET            ATQKEQWFGNRWHEGYRQTPREDSHSTTGTAASAHTSHPMQGRTPSPEDS            SWTDFFPNISHPMGRGHQAGRRMDMDSSHSTTLQPTANPNTGLVEDLDRTGP            LSMTTQQSNSQSFSTSHGLEEDKDHPPTSTLTSSNRNDVTGGRRDPNHSEGS            TTLLEGYTSHPHTKESRTFIPVTSAKTGSFGVTAVTVGDSNSNVNRSLSGDQDT            FHPSGGSHHTHGSSESDGHSHGSQEGGANTTSGPIRTPQIPE</p>
<b>Purity</b>	Greater than 90% as determined by SDS-PAGE.
<b>Source</b>	E.coli
<b>MW</b>	80.2 kDa
<b>Uniprot ID</b>	<b>P16070</b>
<b>Expiration Date</b>	6 months from date of receipt.