

## Product Datasheet

### NALP3/CIAS1 Rabbit Polyclonal Antibody (APC-Cy7) (orb2443833)

|                            |                                                                                  |
|----------------------------|----------------------------------------------------------------------------------|
| <b>Description</b>         | NALP3/CIAS1 Rabbit Polyclonal Antibody (APC-Cy7)                                 |
| <b>Species/Host</b>        | Rabbit                                                                           |
| <b>Reactivity</b>          | Human, Mouse, Rabbit, Rat                                                        |
| <b>Conjugation</b>         | APC/Cy7                                                                          |
| <b>Tested Applications</b> | IF                                                                               |
| <b>Immunogen</b>           | KLH conjugated synthetic peptide derived from human NALP3/CIAS1 (501-600/1036aa) |
| <b>Target</b>              | NLRP3                                                                            |
| <b>Preservatives</b>       | 0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.           |
| <b>Concentration</b>       | 1mg/ml                                                                           |
| <b>Storage</b>             | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.  |
| <b>Note</b>                | For research use only                                                            |
| <b>Isotype</b>             | IgG                                                                              |
| <b>Clonality</b>           | Polyclonal                                                                       |
| <b>Antibody Type</b>       | Primary Antibody                                                                 |
| <b>Dilution Range</b>      | IF=1:100-500                                                                     |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                  |

#### **Biorbyt Ltd.**

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

#### **Biorbyt LLC.**

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)