

Product Datasheet

Human MINPP1 protein (orb244995)

Description	Recombinant human Multiple inositol polyphosphate phosphatase 1
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal GST-tagged
Note	For research use only
Application notes	This is GST-tag protein
Protein Sequence	<p>SLLEPRDPVASSLSPYFGTKTRYEDVNPVLLSGPEAPWRDPELLEGTCTPVQLVA LIRHGTRYPTVKQIRKLRQLHGLLQARGSRDGGASSTGSRDLGAALADWPLWY ADWMDGQLVEKGRQDMRQLALRLASLFPALFSRENYGRLRLITSSKHRCMDSS AAFLQGLWQHYHPGLPPPVDADMEFGPPTVNDKLMRFFDHCEKFLTEVEKNAT ALYHVEAFKTPGEMQNILKVAATLQVPVNDLNADLIQVAFFTCFSDLAIKGVKSP WCDVFDIDDAKVLEYLNDLKQYWKRGYGTINSRSSCTLFQDIFQHLDKAVEQK QRSQPISSPVILQFGHAETLLPLLSLMGYFKDKEPLTAYNYKKQMHRKFRSGLIVP YASNLIFVLYHCENAKTPKEQFRVQMLLNEKVLPLAYSQETVSFYEDLKNHYKDIL QSCQTSEECELARANSTSDDEL</p>
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	79.1 kDa
Uniprot ID	Q9UNW1
Expiration Date	6 months from date of receipt.