

Product Datasheet

Human Btk protein (orb245766)

Description	Recombinant human Tyrosine-protein kinase BTK
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-tagged
Note	For research use only
Application notes	This is His-tag protein
Protein Sequence	AAVILESIFLKRSQKKKTSPLNFKKRLFLTLVHKLSYYEYDFERGRRGSKKGSIDV EKITCVETVVPEKNPPPERQIPRRGEESSEMEQISIIERFPYPFQVVYDEGPLYVFSP TEELRKRWIHLQKNVIRYNSDLVQKYHPCFWIDGQYLCCSQTAKNAMGCQILEN RNGSLKPGSSHRKTKKPLPPTPEEDQILKKPLPPEPAAAPVSTSELKKVVALYDYM PMNANDLQLRKGDEYFILEESNLPWWRRARDKNGQEGYIPSNYVTEAEDSIEMYE WYSKHMTRSQAELLKQEGKEGGFIVRDSSKAGKYTVSVFAKSTGDPQGVIRHY VVCSTPQSQYYLAEKHLFSTIPELINYHQHNSAGLISRLKYPVSQQNKNAPSTAGL GYGSWEIDPKDLTFLKELGTGQFGVVKYGKWARGQYDVAIKMIKEGSMSEDEFIE EAKVMMNLSHEKLVQLYGVCTKQRPFIITEYMANGCLLNLYREMRHRFQTQQLL EMCKDVCEAMEYLESKQFLHRDLAARNCLVNDQGVVKVSDFGLSRYVLDDEYT SSVGSKFPVRWSPPEVLMYSKFSSKSDIWAFGVLMWEIYSLGKMPYERFTNSET AEHIAQGLRLYRPHLASEKVYTIMYSCWHEKADERPTFKILLSNILDVMDEES
Purity	Greater than 90% as determined by SDS-PAGE.
Source	Yeast
MW	78.2 kDa
Uniprot ID	Q06187
Expiration Date	6 months from date of receipt.