

Product Datasheet

E.coli AC1 protein (orb246494)

Description	Recombinant African cassava mosaic virus RepliCation-associated protein
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-SUMO-tagged
Note	For research use only
Application notes	Full length of His-SUMO-tag and expression region is 1-358aa
Protein Sequence	MRTPRFRIQAKNVFLTPKCSIPKEHLLSFIQTLSQLSNPKFIKICRELHQNGEPHL HALIQFEGKITITNNRLFDCVHPSCSTSFHPNIQGAKSSSDVKSYLEDKDGDVIEW GQFQIDGRSARGGQQSANDAYAKALNSGSKSEALNVIRELVPKDFVLQFHNLS NLDRIFQEPPAPYVSPFPCSSFDQVPVEIEEWVADNVDRDAAAPWRPNSIVIEGD SRTGKTIWARSLGPHNYLCGHLDLSPKVFNNAAWYNVIDDVPHYLKHFKEFMG SQRDWQSNTKYGKPVQIKGGIPTIFLCNPGPTSSYKEFLAEEKQEALKAWALKNA IFITLTEPLYSGSNQSHSQTSQEASHPA
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	56.3 kDa
Uniprot ID	P14982
Expiration Date	6 months from date of receipt.