



Product Datasheet

PPP2R1B antibody (orb247581)


Description	Rabbit polyclonal antibody to PPP2R1B
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human PP2A-Abeta/PR65β/PP2A-Aβ/PR65β/PPP2R1B(NP_859050.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	65kDa
Uniprot ID	P30154
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of Jur...



Immunohistochemical staining of mouse lu...



Immunohistochemistry of paraffin-embedde...