



Product Datasheet

RPS14 antibody (orb247595)


| | |
|----------------------------|---|
| Description | Rabbit polyclonal antibody to RPS14 |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | ICC, IF, IHC, WB |
| Immunogen | Recombinant fusion protein of human RPS14(NP_005608.1). |
| Preservatives | PBS with 0.02% sodium azide, 50% glycerol, pH7.3. |
| Concentration | 1 mg/ml |
| Storage | Store at -20°C. Avoid freeze / thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| MW | 16kDa |
| Uniprot ID | P62263 |
| Dilution Range | WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:50-200 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of extracts of var...



Immunofluorescence analysis of MCF7 cell...



Immunohistochemical staining of mouse br...