


Product Datasheet

BCL2 (phospho-S70) antibody (orb248095)

Description	Rabbit polyclonal antibody to BCL2 (phospho-S70)
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	A synthetic phosphorylated peptide around S70 of human Bcl-2(NP_000624.2).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	28kDa
Uniprot ID	P10415
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...