

Product Datasheet

NKIRAS1 Antibody (orb248144)

Description

NKIRAS1 is part of the Ras family and GTPase superfamily

Species/Host

Mouse

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

WB

Applications

Immunogen

A recombinant protein fragment from the central region of NKIRAS1 was used as the immunogen for this antibody.

Preservatives

1X PBS, pH 7.4

Storage

Aliquot the NKIRAS1 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Provided assay concentrations are suggestions only, NKIRAS1 antibody titration may be required for optimal results.

Formula

0.5 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype

Mouse IgG1, kappa

Clonality

Monoclonal

Purity

Protein G purified antibody

Uniprot ID

[Q9NYS0](#)

Hazard

This NKIRAS1 antibody is available for research use only.

Information

Entrez


[28512](#)

Dilution Range

Western blot: 5-10ug/ml

Expiration Date

12 months from date of receipt.



Western blot testing of human samples wi...