

---

## Product Datasheet

### CD90 Antibody / Thy1 (orb248565)

**Description**

CD90 (Thy1) is an 18-35kDa GPI-anchored

many cell types, s...

**Species/Host**

Mouse

**Reactivity**

Human

**Conjugation**

Unconjugated

**Tested Applications**

FACS, IF

**Immunogen**

Purified human brain CD90 / Thy1 was used as the immunogen for this antibody.

**Preservatives**

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

**Storage**

Store the CD90 / Thy1 antibody at 2-8°C.

**Note**

For research use only

**Application notes**

The concentration stated for each application is a general starting point. Variations in protocols, secondaries and substrates may require the CD90 antibody to be titered up or down for optimal performance.

**Formula**

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

**Isotype**

Mouse IgG1, kappa

**Clonality**

Monoclonal

**Purity**

Protein G affinity chromatography

**Hazard Information**

This CD90 / Thy1 antibody is available for research use only.

**Entrez**[7070](#)**Dilution Range**

Flow cytometry: 1-2ug/million cells, Immunofluorescence: 2-4ug/ml

**Expiration Date**

12 months from date of receipt.