

Product Datasheet

Cytomegalovirus Antibody / CMV / PP65 (orb248643)

Description

Ten to forty percent of the patients with

(AIDS) develop cytomega...

Species/Host

Mouse

Reactivity

Virus

Conjugation

Unconjugated

Tested

ELISA, ICC, IF

Applications

Immunogen

Recombinant Cytomegalovirus p65 protein was used as the immunogen for this PP65 antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the PP65 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

The concentration stated for each application is a general starting point. Variations in protocols, secondaries and substrates may require the PP65 antibody to be titrated up or down for optimal performance.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgG1, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Hazard

Information

This PP65 antibody is available for research use only.

Dilution Range

ELISA (order BSA/sodium azide-free format for coating), Immunocytochemistry (Acetone- or paraformaldehyde-fixed): 0.5-1ug/ml for 30 min, Immunofluorescence: 0.5-1ug/ml

Expiration Date

12 months from date of receipt.