

Product Datasheet

PSN632408 (orb251449)

| | |
|------------------------|---|
| Description | PSN632408 is an optimized agonist of GPR119 receptors that shows similar potency to OEA at both reco... |
| Note | For research use only |
| Formula | $C_{18}H_{24}N_4O_4$ |
| MW | 360.41 |
| CAS Number | [857652-30-3] |
| SMILES | <chem>O=C(OC(C)(C)C)N1CCC(CC1)OCC2=NC(C3=CC=NC=C3)=NO2</chem> |
| Expiration Date | 12 months from date of receipt. |