

Product Datasheet

Tandutinib(MLN518) (orb251450)

Description	Tandutinib (MLN518, CT53518) is a potent FLT3 antagonist with IC50 of 0.22 µM, also inhibits PDGFR ...
Note	For research use only
Formula	$C_{31}H_{42}N_6O_4$
MW	562.7
CAS Number	[387867-13-2]
SMILES	<chem>O=C(N1CCN(C2=C(C(C=C3OCCCN4CCCCC4)=NC=N2)C=C3OC)CC1)NC5=CC=C(OC(C)C)C=C5</chem>
Expiration Date	12 months from date of receipt.