

---

## Product Datasheet

### Nrf2 (phospho-Ser40) antibody (orb251473)

Evolution Biocomponents.

Nrf2 (phospho-Ser40) antibody

<b>Description</b>	Nrf2 (phospho-Ser40) antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IHC-P
<b>Immunogen</b>	KLH conjugated synthetic peptide derived from human Nrf2 (phospho-Ser40). Please contact us for the exact immunogen sequence. The peptide is available as orb374857.
<b>Target</b>	Nrf2 (phospho-Ser40)
<b>Form/Appearance</b>	10 mM PBS, 0.02% sodium azide
<b>Concentration</b>	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
<b>Storage</b>	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Purity</b>	Polyclonal antibodies are purified by peptide affinity chromatography
<b>MW</b>	68 kDa
<b>Uniprot ID</b>	<a href="#">Q16236</a> , <a href="#">O54968</a>
<b>NCBI</b>	<a href="#">001138884</a> , <a href="#">21</a> , <a href="#">001145412</a> , <a href="#">11</a>
<b>Entrez</b>	<a href="#">4780</a>
<b>Dilution Range</b>	IHC-P: 1:200
<b>Expiration Date</b>	12 months from date of receipt.