

Product Datasheet

Histone H3(Phospho-Thr32) Rabbit Polyclonal Antibody (orb251700)

Description	Histone H3(Phospho-Thr32) Rabbit Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P68431
Expiration Date	12 months from date of receipt.