

Product Datasheet

Histone H3(Di-methyl-K36) Mouse Monoclonal Antibody (1E6) (orb251795)

Description	Histone H3(Di-methyl-K36) Mouse Monoclonal Antibody (1E6)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Monoclonal
Uniprot ID	P68431
Expiration Date	12 months from date of receipt.