



Product Datasheet

GLP1R Antibody (orb252113)

Description	GLP1R Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P43220
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



Immunohistochemical staining of Human br...