



Product Datasheet

CCL1 Antibody (orb254111)

Description	CCL1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P22362
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of Human colon tissue showing brown staining.



Immunohistochemical analysis of Human colon tissue showing brown staining.