



Product Datasheet

PARP8 Antibody (orb254144)

Description	PARP8 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q8N3A8
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa cell using...



Immunohistochemical analysis of Human CO...