



Product Datasheet

TRIM40 Antibody (orb254223)

Description	TRIM40 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q6P9F5
Expiration Date	12 months from date of receipt.



Western blot analysis of Human normal st...



Immunohistochemical staining of Human li...