



Product Datasheet

VCAM1 antibody (orb254410)

Description	VCAM1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P19320
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...



Immunohistochemistry analysis of paraffi...