



Product Datasheet

LHB antibody (orb254817)


Description	LHB antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P01229
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of K56...



Immunohistochemistry analysis of paraffi...



Immunohistochemistry analysis of paraffi...