

## Product Datasheet

### Duck Hepatitis Virus VP1 Rabbit Polyclonal Antibody (APC-Cy5.5) (orb2550667)

|                             |                                                                                        |
|-----------------------------|----------------------------------------------------------------------------------------|
| <b>Catalog Number</b>       | orb2550667                                                                             |
| <b>Category</b>             | Antibodies                                                                             |
| <b>Description</b>          | Duck Hepatitis Virus VP1 Rabbit Polyclonal Antibody (APC-Cy5.5)                        |
| <b>Clonality</b>            | Polyclonal                                                                             |
| <b>Species/Host</b>         | Rabbit                                                                                 |
| <b>Isotype</b>              | IgG                                                                                    |
| <b>Conjugation</b>          | APC/Cy5.5                                                                              |
| <b>Concentration</b>        | 1mg/ml                                                                                 |
| <b>Buffer/Preservatives</b> | 0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.                 |
| <b>Purification</b>         | Affinity purified by Protein A                                                         |
| <b>Immunogen</b>            | KLH conjugated synthetic peptide derived from Duck Hepatitis Virus VP1 (101-200/240aa) |
| <b>Tested applications</b>  | ICC, IF                                                                                |
| <b>Dilution range</b>       | ICC/IF=1:50-200, IF=1:100-500                                                          |
| <b>Antibody Type</b>        | Primary Antibody                                                                       |
| <b>Storage</b>              | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.        |
| <b>Note</b>                 | For research use only                                                                  |

#### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

#### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)