

Product Datasheet

RFPL4B antibody (Biotin) (orb25589)

Description	RFPL4B antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Ret finger protein-like 4B protein (1-263AA)
Target	RFPL4B
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q6ZWI9
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.