

## Product Datasheet

# Avi Tag Mouse Monoclonal Antibody(6E3) (orb256044)

<b>Description</b>	Avi Tag Mouse Monoclonal Antibody(6E3)
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Storage</b>	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
<b>Note</b>	For research use only
<b>Clonality</b>	Monoclonal
<b>Expiration Date</b>	12 months from date of receipt.