

## Product Datasheet

### Flag-Tag Mouse Monoclonal Antibody (orb256046)

<b>Description</b>	Flag-Tag Mouse Monoclonal Antibody
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IF, IP, WB
<b>Storage</b>	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
<b>Note</b>	For research use only
<b>Clonality</b>	Monoclonal
<b>Expiration Date</b>	12 months from date of receipt.