

Product Datasheet

Porcine Interleukin 17 (IL17) ELISA Kit (orb2560701)

Description	Porcine Interleukin 17 (IL17) ELISA Kit
Reactivity	Porcine
Range	31.2-2000pg/mL
Storage	Store at 4°C & -20°C for 12 months.
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sensitivity	12.5pg/mL
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)