


Product Datasheet

AmCyan Mouse Polyclonal Antibody (orb256071)

Description	AmCyan Mouse Polyclonal Antibody
Species/Host	Mouse
Reactivity	Canine, Gallus, Hamster, Human, Monkey, Mouse, Porcine, Rabbit, Rat
Conjugation	Unconjugated
Tested Applications	WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.



Western blot analysis of AmCyan recombin...