


## Product Datasheet

### V5 tag antibody (orb256432)

<b>Description</b>	Mouse monoclonal antibody to V5 tag
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IF, IHC-Fr, IP, RIA, WB
<b>Immunogen</b>	Pk (V5) Epitope Tag (GKPIPPLLGLDST).
<b>Target</b>	V5 epitope tag
<b>Preservatives</b>	PBS with 0.02% Proclin 300.
<b>Concentration</b>	1 mg/ml
<b>Storage</b>	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Purity</b>	Purified
<b>Expiration Date</b>	12 months from date of receipt.



ELISA of anti-epitope tag antibodies on

...