

Product Datasheet

Mouse IgG2B kappa isotype control (orb256459)

Description	Mouse IgG2B kappa isotype control
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	Control
Target	Unknown Specificity (Isotype control)
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.