

## Product Datasheet

### Rat IgG1 kappa isotype control (orb256462)

<b>Description</b>	Rat IgG1 kappa isotype control
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	Control
<b>Target</b>	Unknown Specificity (Isotype control)
<b>Preservatives</b>	PBS with 0.02% Proclin 300.
<b>Concentration</b>	1 mg/ml
<b>Storage</b>	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Purity</b>	Purified
<b>Expiration Date</b>	12 months from date of receipt.