

Product Datasheet

E3 Ligase Ligand-linker Conjugate 91 (orb2566355)

Description	E3 Ligase Ligand-linker Conjugate 91, a conjugate of an E3 ligase ligand and linker, consists of (S,R,S)-AHPC and the corresponding linker. (S,R,S)-AHPC is a VH032-based von Hippel-Lindau (VHL) ligand that recruits VHL proteins. This conjugate serves as a key intermediate for synthesizing complete PROTAC molecules.
Storage	-20°C
Note	For research use only
Formula	C ₃₅ H ₅₂ N ₆ O ₇ S
Purity	98.00%
MW	700.89
CAS Number	2229713-40-8
SMILES	<chem>C([C@@H])(NC(COCCN1CCN(C(OC(C)(C)C)=O)CC1)=O)[C@@](C)(C)C(=O)N2[C@H](C(NCC3=CC=C(C=C3)C4=C(C)N=CS4)=O)C[C@@H](O)C2</chem>
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)