

Product Datasheet

E3 Ligase Ligand-linker Conjugate 5 (orb2566361)

Description	E3 Ligase Ligand-linker Conjugate 5, a combination of an E3 ligase ligand and linker, comprises Thalidomide and its specific Linker. This conjugate acts as a Cereblon ligand to attract the CRBN protein and functions as a crucial intermediate in the production of full PROTAC molecules.
Storage	-20°C
Note	For research use only
Formula	$C_{29}H_{37}N_5O_6$
Purity	98.00%
MW	551.63
CAS Number	2839670-45-8
SMILES	<chem>C(OC(C)(C)C)(=O)N1CC2(C1)CN(CC3CCN(CC3)C=4C=C5C(=CC4)C(=O)N(C5=O)C6C(=O)NC(=O)CC6)C2</chem>
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)