

Product Datasheet

E3 Ligase Ligand-linker Conjugate 45 (orb2566363)

Description	E3 Ligase Ligand-linker Conjugate 45, a conjugate of E3 ligase ligand and linker, consists of Thalidomide and the corresponding Linker. It can serve as a Cereblon ligand to recruit CRBN protein and acts as a key intermediate for synthesizing complete PROTAC molecules.
Storage	-20°C
Note	For research use only
Formula	C ₂₅ H ₃₁ N ₅ O ₆
Purity	98.00%
MW	497.54
CAS Number	2871775-00-5
SMILES	<chem>O=C1N(C(=O)C=2C1=CC(=CC2)N3CCN(CC4CCN(CC(O)=O)CC4)CC3)C5C(=O)NC(=O)CC5</chem>
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)