

Product Datasheet

Ac-WEHD-pNA (orb2567076)

Description	Ac-WEHD-pNA serves as a colorimetric substrate for caspases 1, 4, 5, and 14, widely used to assess caspase activity [1].
Storage	-20°C
Note	For research use only
Formula	C ₃₄ H ₃₇ N ₉ O ₁₁
Purity	98.00%
MW	747.71
CAS Number	1135248-58-6
SMILES	<chem>C([C@@H](C(N[C@H](C(N[C@@H](CC1=CN=CN1)C(N[C@H](C(NC2=CC=C(N(=O)C=C2)=O)C=C2)=O)CC(O)=O)=O)=O)CCC(O)=O)=O)NC(C)=O)C=3C=4C(NC3)=CC=CC4</chem>
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)