

Product Datasheet

beta-catenin (phospho-S33) peptide (orb256772)

Description	The peptide is used to anti beta-catenin (phospho-S33) antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rabbit, Rat, Sheep
Tested Applications	Blocking
Target	CTNNB1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb256482 (beta-catenin (phospho-S33) antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9WU82 , Q02248 , P35222
Entrez	12387 , 84353 , 1499
Expiration Date	12 months from date of receipt.